

Solar Cell with Motor


Assembly Instructions

Listed below are the components of the Solar Cell and instructions for preparing a globe for the seasons investigation.

Parts Checklist

The following items are provided with the Solar Cell & Motor:

- wires with banana plugs (red & black)
- alligator clips (2)
- round Velcro dots (8)
- Solar Cell
- motor
- propeller



Attaching Velcro dots to a Globe (for Seasons Investigation)

1. Place this piece of paper on the globe, matching the Meridian to the dotted line below.
(The dotted line aligns with the Meridian. The arrows indicate the placement of Velcro dots.)
2. Affix one dot on the globe, at each arrow.
(Place 2 dots at equator, 2 dots at 45° North, 2 dots at 45° South.)

For detailed instructions on these investigations, consult your CPO Science Teacher's Guide.

For technical assistance, please call 866.588.6951.

